

Procedure to start the run

- 1) Clear the spectra server in Vigru
- 2) Ensure that online, the spectra spy, is started. To ensure that, from a terminal:
 - > *Online*
 - > *start_online*
- 3) Start the run from Run Control (click **start**)
- 4) Check that the scalers are counting properly
- 5) Take a screen-shot of the scalers. Name it: Scalers_start_runNbr.png
- 6) Write in the elog the relevant info about the run:
 - User: **e748**
 - Password: **Belise!**

Author: The persons making the entry in the elog

ExpNbr: e744

Type: Run (unless particular type of run)

Category: General

Status: Ongoing

Target-Source: CH2, Carbon, etc

Subject: Run number and target nr (from 1 -top to 6-bottom)

entry text body:

Run number

Beam used (14N, 14O)

start time

stop time

Target used in the run

Persons in shift

Add the file: Scalers_start_runNbr.png

Example:

Author: Iulian STEFAN

ExpNbr: e744

Type: Run

Category: General

Status: Ongoing

Target-Source: CH2

Subject: Run 500 target nr. 2 (CH2 100 um)

entry text body:

Run 500

Beam 14O

start time 10h30m10s 26/05/2018

stop time

Target nr 2 (CH2 100 um)

Iulian STEFAN, Megumi NIIKURA, Fairouz HAMACHE

added file Scalers_start_run500.png